

Comparison of lens orientation stability of two daily disposable silicone hydrogel toric lenses

Jill Woods,

Amir Moezzi, Jalaiah Varikooty, Sarah Guthrie, Lyndon Jones



UNIVERSITY OF WATERLOO
FACULTY OF SCIENCE
School of Optometry & Vision Science



*Disclosure: study sponsored by CooperVision
CORE: funded by all major CL manufacturers*



Purpose & background

- There are several daily disposable (DD) silicone hydrogel (SH) options to correct astigmatism
- Soft toric lenses require stabilization designs to maintain the correct orientation of the cylindrical power.
- This study compared orientation stability of 2 x DDSHs of different stabilization designs.


MDT: MyDay toric

(stenfilcon A; Optimised Toric Lens Geometry™, CooperVision)

AO1A: Acuvue® Oasys 1-Day with Hydraluxe™ for Astigmatism

(senofilcon A; Blink Stabilized® design, Johnson & Johnson)






MDT



- 1. Uniform horizontal ISO thickness**
 - Improves lens stability and reduces rotation for better visual acuity*
- 2. Smooth, continuous ballast**
 - Maximizes comfort* during eyelid interaction
- 3. Larger toric optic zone**
 - Provides clear visual acuity*
- 4. Optimized ballast toric design**
 - Provides a stable, comfortable fit*

A01A



BLINK STABILIZED™

Balanced with four
Zones of stability





Methods

- 1 site - CORE, 20 subjects, habitual soft toric wearers
- 1-week cross-over, randomized lens order, double masked
- Key outcomes, comparing data from 1-week visits :
 - LogMAR acuity & lens orientation with
 - Head straight
 - Head tilted 30° to right
 - Head tilted 30° to left
 - Rotational recovery (# degrees lens rotates back in 10 blinks after 30° forced lens rotation)






Methods

HEAD TILT	ROTATIONAL RECOVERY
	
<ul style="list-style-type: none"> • Tilt head 30°, wait for >1 minute. • Measure lens orientation relative to 'new' vertical position (β). • If $\beta = 18^\circ$, lens rotates 12° from habitual. • Best case is $\beta = 30^\circ$ ie. no rotation on head tilt 	<ul style="list-style-type: none"> • Rotate lens 30° from habitual position • After 10 normal blinks, measure degrees recovered. • Best recovery = 30°, returned to habitual position • Worst recovery = 0°



Results: RE mean data

			
LOGMAR ACUITY	MDT -0.07 AO1A -0.08 p=0.45 ½ letter diff	MDT -0.03 AO1A -0.04 p=0.07 ½ letter diff	MDT -0.03 AO1A -0.06 p=0.049 1½ letter diff
LENS ORIENTATION relative to vertical	MDT 6° AO1A 3° p=0.01 (Absolute values)	β: MDT 18° AO1A 17° p=0.41 Lens rotates temporal by MDT 12° AO1A 13°	β : MDT 17° AO1A 19° p=0.20 Lens rotates nasal by MDT 13° AO1A 11°
ROTATION RECOVERY Degrees rotated back after 10 blinks (full recovery = 30°)	Nasal: MDT 27° AO1A 19° p<0.01 Temp: MDT 22° AO1A 17° p=0.01	n/a	n/a



Conclusions

MDT and AO1A have DIFFERENT STABILISATION designs:

- BOTH lenses performed well for visual acuity & lens orientation
- BOTH lenses rotated similar amount to follow the head tilt, <50% of head tilt angle
- MDT showed a faster return when displaced away from habitual position

